



Anesthesia accessory

CATALOGUE

2018.08

www.mindray.com

P/N:ENG-Anesthesia accessory catalogue-210210X58P-20180821
©2018 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.



Anesthesia Mask



Disposable Anesthesia Mask

- Single patient use to avoid cross-infection
- Ergonomic design for convenient one-hand operation
- Adjustable air cushion design for optimal comfort and tightness
- Transparent design for optimal visibility
- Materials in accordance with ISO-10993:10993 biological compatibility requirements
- Complete range of models with color-coded hook for easier identification

Picture	Part No.	Description	Apply to
	040-001817-00	Aircushion mask, disposable, size #0, infant	All
	040-001818-00	Aircushion mask, disposable, size #1, infant large	All
	040-001819-00	Aircushion mask, disposable, size #2, child	All
	040-001820-00	Aircushion mask, disposable, size #3, child large	All

Picture	Part No.	Description	Apply to
	040-001821-00	Aircushion mask, disposable, size #4, adult	All

	040-001822-00	Aircushion mask, disposable, size #5, adult large	All
---	---------------	---	-----

Reusable Anesthesia Mask

- Ergonomic design for convenient one-hand operation
- Transparent design for optimal visibility
- Materials in accordance with ISO-10993:10993 biological compatibility requirements
- Full range of models for different clinical applications



Picture	Part No.	Description	Apply to
	040-001835-00	Silicon face mask, reusable, size #0, infant	All

	040-001836-00	Silicon face mask, reusable, size #1, infant large	All
---	---------------	--	-----

	040-001837-00	Silicon face mask, reusable, size #2, child	All
--	---------------	---	-----

	040-001841-00	Silicon face mask, reusable, size #3, child large	All
---	---------------	---	-----


Vaporizer



- Large capacity up to 360 ml with dry wick
- Highly compatible supporting up to 3 different anesthetic agents
- Automatic flow, temperature, pressure compensation
- Maintenance-free design

Mindray V60 Vaporizer

Picture	Part No.	Description	Apply to
	115-005348-00	Mindray V60 vaporizer, Isoflurane, Pour-fill	All
	115-005345-00	Mindray V60 vaporizer, Isoflurane, Key-fill, with Key-fill adaptor	All
	115-005349-00	Mindray V60 vaporizer, Sevoflurane, Pour-fill	All
	115-005346-00	Mindray V60 vaporizer, Sevoflurane, Key-fill, with Key-fill adaptor	All

Picture	Part No.	Description	Apply to
	040-000065-00	Key-fill Adaptor for filling the vaporizer, Isoflurane	All
	040-000066-00	Key-fill Adaptor for filling the vaporizer, Sevoflurane	All
	115-026747-00	Quik-fill Adaptor for filling the vaporizer, Sevoflurane	All
	040-000067-00	Quik-fill Drain Funnel Adaptor for draining the vaporizer, Sevoflurane	All

Flow sensor

Picture	Part No.	Description	Apply to
	0601-30-69700	Inspiration flow sensor, 1 pcs	WATO EX-65 Pro/ EX-55 Pro/ EX-65/EX-55/ EX-35/EX-30/ EX-20
	0601-30-78894	Expiration flow sensor, 1 pcs	
	115-001366-00	Flow sensor kit, including: - Inspiration flow sensor, 1 pcs 0601-30-69700 - Expiration flow sensor, 1 pcs 0601-30-78894	
	115-008262-00	Inspiration Flow sensor, 1 pcs	A7/A5
	115-008263-00	Expiration Flow sensor, 1 pcs	
	115-008264-00	Flow sensor kit, including: - Inspiration Flow sensor, 1 pcs 115-008262-00 - Expiration Flow sensor, 1 pcs 115-008263-00	
	115-041507-00	Inspiration flow sensor, autoclavable, 1 pcs	A7/A5/ WATO EX-65 Pro/ EX-55 Pro/ EX-65/ EX-55/ EX-35
	115-041508-00	Expiration flow sensor, autoclavable, 1 pcs	
	115-041519-00	Flow sensor kit, autoclavable, including: - Inspiration flow sensor, 1 pcs 115-041507-00 - Expiration flow sensor, 1 pcs 115-041508-00	

